


SENSABLACK **FLAT**

FLAT SURFACE pH ELECTRODE

INDIA'S FIRST ELECTROCHEMISTRY SENSOR WITH:

- Ergonomic "2 tone Cap" design
- Unbreakable Electrode
- H Gel + Filled



TECHNICAL SPECIFICATIONS:

pH Range	0 - 14
Operating Temperature	-5 to 100°C
Slope	59±3 mV (95-105%) at 25°C
Zero Point	pH 7
Asymmetry Potential	0±25 mV
Membrane Glass	Green
Membrane Glass Impedance	350 - 500 MΩ
Diaphragm	Teflon
Junctions	Single
Reference System	Ag/ AgCl
Reference Electrolyte	H Gel +
pH Storage Solution	Saturated KCl/ TOSHCON storage solution
Shape of Membrane	Flat
Diameter	12 mm
Shaft Length	120 mm
Shaft Material	Polycarbonate
Connector type	Cap/ TRC* with BNC with Low-Noise Coaxial Cable or S8 or any other connector required.
Suitability	Suitable/ compatible with mostly all make instruments.

ORDERING CODE

Cap / Connector

Cap	-	2 tone
S8	-	S8

Sensablack Flat-X-X-X-X

Junction:

Single	-	-
--------	---	---

Sensablack Flat-X-X-X-X

Temp. Compensator:

Without	-	-
Pt100	-	T
Pt1000	-	1T

Sensablack Flat-X-X-X-X

Cable Length for Cap Electrode:

0.75	-	-
1 Mtrs.	-	1
2 Mtrs.	-	2
3 Mtrs.	-	3

Sensablack Flat-X-X-X-X

OTHER INDUSTRY:

Pharma, Petrochemical, Electroplating, Food, Chemical, Fertilizer, Dairy, Oil, Drinking Water, Ground Water, Waste Water.

Manufactured by:

TOSHNIWAL INSTRUMENTS MANUFACTURING PVT. LTD.,
P.O. Gagwana - 305023, Distt. Ajmer (Raj.) INDIA, PHONE: +91-145-2971131/2//3
CUSTOMER CARE: +91-9829112345

EMAIL: info@toshcon.com, support@toshcon.com

REGD OFFICE: 401, Manish Chamber, Sonawala Cross Lane, Goregaon (E) Mumbai-400 063

WEBSITE: www.toshcon.com